

Technical Data Muffle Furnaces up to 1100°C

Muffle Furnaces KL 03/11

Max. temperature: 1100°C
Inside dimensions (mm): W 190 x D 140 x H 100
Volume: 2,7 l
Outside dimensions (mm): W 437 x D 444 x H 489
Power: 1,2 kW
Voltage: 230 1/N V
Weight: 20 kg

Muffle Furnaces KL 05/11

Max. temperature: 1100°C
Inside dimensions (mm): W 240 x D 170 x H 130
Volume: 5,3 l
Outside dimensions (mm): W 487 x D 462 x H 519
Power: 2,4 kW
Voltage: 230 1/N V
Weight: 35 kg

Muffle Furnaces KL 09/11

Max. temperature: 1100°C
Inside dimensions (mm): W 240 x D 260 x H 170
Volume: 10,6 l
Outside dimensions (mm): W 492 x D 537 x H 558
Power: 3,0 kW
Voltage: 230 1/N V
Weight: 45 kg

Muffle Furnaces KL 15/11

Max. temperature: 1100°C
Inside dimensions (mm): W 250 x D 370 x H 170
Volume: 15,7 l
Outside dimensions (mm): W 499 x D 648 x H 558
Power: 3,5 kW
Voltage: 230 1/N V
Weight: 50 kg