

Technical Data Split-Type Tube Furnaces up to 1100°C

Split-type Tube Furnaces ROK 70/250/11

Max. temperature:	1100°C
Heated length:	250 mm
Temp. const. zone:	80 mm
Outer dimensions (mm):	W 400 x D 500 x H 510
Max. Ø pipe intake:	≤ 70 mm
Power:	1,6 kW
Voltage:	230 1/N V

Split-type Tube Furnaces ROK 70/500/11

Max. temperature:	1100°C
Heated length:	500 mm
Temp. const. zone:	160 mm
Outer dimensions (mm):	W 650 x D 500 x H 510
Max. Ø pipe intake:	≤ 70 mm
Power:	2,4 kW
Voltage:	230 1/N V

Split-type Tube Furnaces ROK 100/250/11

Max. temperature:	1100°C
Heated length:	250 mm
Temp. const. zone:	80 mm
Outer dimensions (mm):	W 400 x D 500 x H 510
Max. Ø pipe intake:	≤ 100 mm
Power:	2,4 kW
Voltage:	230 1/N V

Split-type Tube Furnaces ROK 100/500/11

Max. temperature:	1100°C
Heated length:	500 mm
Temp. const. zone:	160 mm
Outer dimensions (mm):	W 650 x D 500 x H 510
Max. Ø pipe intake:	≤ 100 mm
Power:	3,0 kW
Voltage:	230 1/N V

Split-type Tube Furnaces ROK 150/250/11

Max. temperature:	1100°C
Heated length:	250 mm
Temp. const. zone:	80 mm
Outer dimensions (mm):	W 400 x D 500 x H 510
Max. Ø pipe intake:	≤ 130 mm
Power:	3,0 kW
Voltage:	230 1/N V

Split-type Tube Furnaces ROK 150/500/11

Max. temperature:	1100°C
Heated length:	500 mm
Temp. const. zone:	160 mm
Outer dimensions (mm):	W 650 x D 500 x H 510
Max. Ø pipe intake:	≤ 130 mm
Power:	3,5 kW
Voltage:	230 1/N V

Split-type Tube Furnaces ROK 200/250/11

Max. temperature:	1100°C
Heated length:	250 mm
Temp. const. zone:	80 mm
Outer dimensions (mm):	W 400 x D 500 x H 510
Max. Ø pipe intake:	≤ 160 mm
Power:	3,0 kW
Voltage:	230 1/N V

Split-type Tube Furnaces ROK 200/500/11

Max. temperature:	1100°C
Heated length:	500 mm
Temp. const. zone:	160 mm
Outer dimensions (mm):	W 650 x D 500 x H 510
Max. Ø pipe intake:	≤ 160 mm
Power:	3,5 kW
Voltage:	230 1/N V

Split-type Tube Furnaces ROK 250/400/11

Max. temperature:	1100°C
Heated length:	400 mm
Temp. const. zone:	130 mm
Outer dimensions (mm):	W 555 x D 740 x H 710
Max. Ø pipe intake:	≤ 200 mm
Power:	6,0 kW
Voltage:	400 3/N V

Split-type Tube Furnaces ROK 300/400/11

Max. temperature:	1100°C
Heated length:	400 mm
Temp. const. zone:	130 mm
Outer dimensions (mm):	W 555 x D 740 x H 710
Max. Ø pipe intake:	≤ 250 mm
Power:	9,0 kW
Voltage:	230 1/N V