

## Technical Data Tube Furnaces up to 1800°C

### Tube Furnaces ROHT 40/200/18

<b>Max. temperature:</b>	<b>1800°C</b>
<b>Heated length:</b>	<b>200 mm</b>
<b>Outer dimensions (mm):</b>	<b>W 460 x D 530 x H 830</b>
<b>Tube- Ø:</b>	<b>40 mm</b>
<b>Power:</b>	<b>3,0 kW</b>
<b>Voltage:</b>	<b>230 1/N V</b>

### Tube Furnaces ROHT 40/300/18

<b>Max. temperature:</b>	<b>1800°C</b>
<b>Heated length:</b>	<b>300 mm</b>
<b>Outer dimensions (mm):</b>	<b>W 560 x D 530 x H 830</b>
<b>Tube- Ø:</b>	<b>40 mm</b>
<b>Power:</b>	<b>4,0 kW</b>
<b>Voltage:</b>	<b>400 3/N V</b>

### Tube Furnaces ROHT 40/400/18

<b>Max. temperature:</b>	<b>1800°C</b>
<b>Heated length:</b>	<b>400 mm</b>
<b>Outer dimensions (mm):</b>	<b>W 660 x D 530 x H 830</b>
<b>Tube- Ø:</b>	<b>40 mm</b>
<b>Power:</b>	<b>4,8 kW</b>
<b>Voltage:</b>	<b>400 3/N V</b>

### Tube Furnaces ROHT 40/600/18

<b>Max. temperature:</b>	<b>1800°C</b>
<b>Heated length:</b>	<b>600 mm</b>
<b>Outer dimensions (mm):</b>	<b>W 860 x D 530 x H 830</b>
<b>Tube- Ø:</b>	<b>40 mm</b>
<b>Power:</b>	<b>8,6 kW</b>
<b>Voltage:</b>	<b>400 3/N V</b>

#### Tube Furnaces ROHT 50/200/18

**Max. temperature:** 1800°C  
**Heated length:** 200 mm  
**Outer dimensions (mm):** W 460 x D 530 x H 830  
**Tube- Ø:** 50 mm  
**Power:** 3,0 kW  
**Voltage:** 230 1/N V

#### Tube Furnaces ROHT 50/300/18

**Max. temperature:** 1800°C  
**Heated length:** 300 mm  
**Outer dimensions (mm):** W 560 x D 530 x H 830  
**Tube- Ø:** 50 mm  
**Power:** 4,0 kW  
**Voltage:** 400 3/N V

#### Tube Furnaces ROHT 50/400/18

**Max. temperature:** 1800°C  
**Heated length:** 400 mm  
**Outer dimensions (mm):** W 660 x D 530 x H 830  
**Tube- Ø:** 50 mm  
**Power:** 4,8 kW  
**Voltage:** 400 3/N V

#### Tube Furnaces ROHT 50/600/18

**Max. temperature:** 1800°C  
**Heated length:** 600 mm  
**Outer dimensions (mm):** W 860 x D 530 x H 830  
**Tube- Ø:** 50 mm  
**Power:** 8,6 kW  
**Voltage:** 400 3/N V

#### Tube Furnaces ROHT 75/200/18

**Max. temperature:** 1800°C  
**Heated length:** 200 mm  
**Outer dimensions (mm):** W 460 x D 530 x H 830  
**Tube- Ø:** 75 mm  
**Power:** 3,0 kW  
**Voltage:** 230 1/N V

#### Tube Furnaces ROHT 75/300/18

**Max. temperature:** 1800°C  
**Heated length:** 300 mm  
**Outer dimensions (mm):** W 560 x D 530 x H 830  
**Tube- Ø:** 75 mm  
**Power:** 4,0 kW  
**Voltage:** 400 3/N V

#### Tube Furnaces ROHT 75/400/18

**Max. temperature:** 1800°C  
**Heated length:** 400 mm  
**Outer dimensions (mm):** W 660 x D 530 x H 830  
**Tube- Ø:** 75 mm  
**Power:** 4,8 kW  
**Voltage:** 400 3/N V

#### Tube Furnaces ROHT 75/600/18

**Max. temperature:** 1800°C  
**Heated length:** 600 mm  
**Outer dimensions (mm):** W 860 x D 530 x H 830  
**Tube- Ø:** 75 mm  
**Power:** 8,6 kW  
**Voltage:** 400 3/N V