

Technical Data Chamber Furnaces (closed Ceramic Muffle) up to 1300°C

Chamber Furnaces KLS 03/10/M

Max. temperature:	1000°C
Inside dimensions (mm):	W 140 x D 200 x H 110
Volume:	3 l
Outside dimensions (mm):	W 320 x D 350 x H 410
Power:	1,5 kW
Voltage:	230 1/N V
Weight:	20 kg

Chamber Furnaces KLS 07/11/M

Max. temperature:	1150°C
Inside dimensions (mm):	W 210 x D 280 x H 110
Volume:	7 l
Outside dimensions (mm):	W 550 x D 580 x H 650
Power:	2,7 kW
Voltage:	230 1/N V