

Technical Data High Temperature Furnaces / Bench Top up to 1750°C

Hightemperature Furnaces HTL 02/17

Max. temperature: 1750°C
Inside dimensions (mm): W 110 x D 120 x H 150
Volume: 2 l
Outside dimensions (mm): W 740 x D 440 x H 630
Power: 1,5 kW
Voltage: 230 1/N V
Weight: 76 kg

Hightemperature Furnaces HTL 04/17

Max. temperature: 1750°C
Inside dimensions (mm): W 150 x D 150 x H 150
Volume: 3,5 l
Outside dimensions (mm): W 740 x D 440 x H 630
Power: 3,0 kW
Voltage: 230 1/N V
Weight: 83 kg

Hightemperature Furnaces HTL 10/17

Max. temperature: 1750°C
Inside dimensions (mm): W 200 x D 250 x H 200
Volume: 10 l
Outside dimensions (mm): W 790 x D 540 x H 680
Power: 4,0 kW
Voltage: 400 3/N V
Weight: 97 kg

Hightemperature Furnaces HTL 16/17

Max. temperature: 1750°C
Inside dimensions (mm): W 200 x D 300 x H 260
Volume: 15,6 l
Outside dimensions (mm): W 830 x D 550 x H 735
Power: 8 kW
Voltage: 400 3/N V
Weight: 175 kg

Hightemperature Furnaces HTL 20/17

Max. temperature: 1750°C
Inside dimensions (mm): W 250 x D 320 x H 260
Volume: 21 l
Outside dimensions (mm): W 880 x D 570 x H 735
Power: 8 kW
Voltage: 400 3/N V
Weight: 210 kg