

Technical Data Vacuum Drying Ovens

Vacuum Drying Ovens KTL 20/02/V

Max. temperature:	200°C
Inside dimensions (mm):	W 340 x D 260x H 300
Volume:	26,0 l
Outside dimensions (mm):	W 560 x D 490x H 700
Power:	1,0 kW
Voltage:	230 1/N V
Weight:	65 kg

Vacuum Drying Ovens KTL 60/02/V

Max. temperature:	200°C
Inside dimensions (mm):	W 400 x D 320 x H 430
Volume:	55,0 l
Outside dimensions (mm):	W 620 x D 550 x H 830
Power:	1,2 kW
Voltage:	230 1/N V
Weight:	98 kg

Vacuum Drying Ovens KTL 120/02/V

Max. temperature:	200°C
Inside dimensions (mm):	W 540 x D 410 x H 480
Volume:	106,0 l
Outside dimensions (mm):	W 760 x D 640 x H 880
Power:	1,8 kW
Voltage:	230 1/N V
Weight:	130 kg