

Technical Data High-Temperature Drying Ovens

High-Temperature Drying Ovens KU 40/04/A

Max. temperature:	450°C
Inside dimensions (mm):	W 300 x D 400 x H 300
Volume:	35,0 l
Outside dimensions (mm):	W 570 x D 1020 x H 1305
Power:	3,2 kW
Voltage:	230 1/N V
Weight:	115 kg

High-Temperature Drying Ovens KU 70/04/A

Max. temperature:	450°C
Inside dimensions (mm):	W 350 x D 500 x H 400
Volume:	70,0 l
Outside dimensions (mm):	W 620 x D 1120 x H 1405
Power:	4,3 kW
Voltage:	400 3/N V
Weight:	130 kg

High-Temperature Drying Ovens KU 140/04/A

Max. temperature:	450°C
Inside dimensions (mm):	W 450 x D 600 x H 500
Volume:	140,0 l
Outside dimensions (mm):	W 720 x D 1220 x H 1505
Power:	5,8 kW
Voltage:	400 3/N V
Weight:	205 kg

High-Temperature Drying Ovens KU 270/04/A

Max. temperature:	450°C
Inside dimensions (mm):	W 600 x D 750 x H 600
Volume:	270,0 l
Outside dimensions (mm):	W 870 x D 1390 x H 1605
Power:	9,8 kW
Voltage:	400 3/N V
Weight:	370 kg

High-Temperature Drying Ovens KU 540/04/A

Max. temperature:	450°C
Inside dimensions (mm):	W 750 x D 900 x H 800
Volume:	540,0 l
Outside dimensions (mm):	W 1020 x D 1540 x H 1805
Power:	13,8 kW
Voltage:	400 3/N V
Weight:	540 kg