

Technical Data High-Temperature Drying Ovens

High-Temperature Drying Ovens KU 15/06/A

Max. temperature:	650°C
Inside dimensions (mm):	W 300 x D 350 x H 150
Volume:	15,0 l
Outside dimensions (mm):	W 500 x D 900 x H 440
Power:	2,6 kW
Voltage:	230 1/N V
Weight:	63 kg

High-Temperature Drying Ovens KU 40/06/A

Max. temperature:	650°C
Inside dimensions (mm):	W 300 x D 400 x H 300
Volume:	35,0 l
Outside dimensions (mm):	W 540 x D 1050 x H 1270
Power:	4,3 kW
Voltage:	400 3/N V
Weight:	125 kg

High-Temperature Drying Ovens KU 70/06/A

Max. temperature:	650°C
Inside dimensions (mm):	W 350 x D 500 x H 400
Volume:	70,0 l
Outside dimensions (mm):	W 590 x D 1150 x H 1370
Power:	8,3 kW
Voltage:	400 3/N V
Weight:	140 kg

High-Temperature Drying Ovens KU 140/06/A

Max. temperature:	650°C
Inside dimensions (mm):	W 450 x D 600 x H 500
Volume:	135,0 l
Outside dimensions (mm):	W 690 x D 1250 x H 1470
Power:	12,3 kW
Voltage:	400 3/N V
Weight:	220 kg

High-Temperature Drying Ovens KU 270/06/A

Max. temperature:	650°C
Inside dimensions (mm):	W 600 x D 750 x H 600
Volume:	270,0 l
Outside dimensions (mm):	W 840 x D 1450 x H 1620
Power:	16,8 kW
Voltage:	400 3/N V
Weight:	480 kg

High-Temperature Drying Ovens KU 540/06/A

Max. temperature:	650°C
Inside dimensions (mm):	W 750 x D 900 x H 800
Volume:	540,0 l
Outside dimensions (mm):	W 990 x D 1600 x H 1820
Power:	24,8 kW
Voltage:	400 3/N V
Weight:	560 kg