

Technical Data High-Temperature Drying Ovens

High-Temperature Drying Ovens KU 40/07/A

Max. temperature:	750°C
Inside dimensions (mm):	W 300 x D 400 x H 300
Volume:	35,0 l
Outside dimensions (mm):	W 610 x D 1050 x H 1410
Power:	5,5 kW
Voltage:	400 3/N V
Weight:	140 kg

High-Temperature Drying Ovens KU 70/07/A

Max. temperature:	750°C
Inside dimensions (mm):	W 350 x D 500 x H 400
Volume:	70,0 l
Outside dimensions (mm):	W 660 x D 1150 x H 1510
Power:	10,7 kW
Voltage:	400 3/N V
Weight:	170 kg

High-Temperature Drying Ovens KU 140/07/A

Max. temperature:	750°C
Inside dimensions (mm):	W 450 x D 600 x H 500
Volume:	135,0 l
Outside dimensions (mm):	W 760 x D 1260 x H 1610
Power:	14,3 kW
Voltage:	400 3/N V
Weight:	290 kg

High-Temperature Drying Ovens KU 270/07/A

Max. temperature:	750°C
Inside dimensions (mm):	W 600 x D 750 x H 600
Volume:	270,0 l
Outside dimensions (mm):	W 950 x D 1460 x H 1750
Power:	21,7 kW
Voltage:	400 3/N V
Weight:	590 kg

High-Temperature Drying Ovens KU 540/07/A

Max. temperature:	750°C
Inside dimensions (mm):	W 750 x D 900 x H 800
Volume:	540,0 l
Outside dimensions (mm):	W 1100 x D 1610 x H 1950
Power:	28,7 kW
Voltage:	400 3/N V
Weight:	700 kg