

Technical Data High Temperature Chamber Furnaces

High Temperature Chamber Furnaces HKT 16/18

Max. temperature:	1800°C
Inside dimensions (mm):	W 200 x D 300 x H 260
Volume:	16 l
Outside dimensions (mm):	W 820 x D 655 x H 1650
Power:	11,0 kW
Voltage:	400 3/N V
Weight:	310 kg

High Temperature Chamber Furnaces HKT 20/18

Max. temperature:	1800°C
Inside dimensions (mm):	W 250 x D 320 x H 260
Volume:	20 l
Outside dimensions (mm):	W 870 x D 675 x H 1650
Power:	11,0 kW
Voltage:	400 3/N V
Weight:	325 kg

High Temperature Chamber Furnaces HKT 40/18

Max. temperature:	1800°C
Inside dimensions (mm):	W 300 x D 350 x H 350
Volume:	40 l
Outside dimensions (mm):	W 920 x D 705 x H 1740
Power:	11,0 kW
Voltage:	400 3/N V
Weight:	370 kg

High Temperature Chamber Furnaces HKT 50/18

Max. temperature:	1800°C
Inside dimensions (mm):	W 250 x D 550 x H 350
Volume:	50 l
Outside dimensions (mm):	W 870 x D 905 x H 1740
Power:	18,0 kW
Voltage:	400 3/N V
Weight:	455 kg

High Temperature Chamber Furnaces HKT 70/18

Max. temperature:	1800°C
Inside dimensions (mm):	W 400 x D 400 x H 400
Volume:	65 l
Outside dimensions (mm):	W 1100 x D 835 x H 1740
Power:	16,0 kW
Voltage:	400 3/N V
Weight:	540 kg

High Temperature Chamber Furnaces HKT 100/18

Max. temperature: 1800°C
Inside dimensions (mm): W 400 x D 600 x H 400
Volume: 100 l
Outside dimensions (mm): W 1100 x D 1035 x H 1790
Power: 21,0 kW
Voltage: 400 3/N V
Weight: 620 kg

High Temperature Chamber Furnaces HKT 130/18

Max. temperature: 1800°C
Inside dimensions (mm): W 400 x D 800 x H 400
Volume: 130 l
Outside dimensions (mm): W 1100 x D 1235 x H 1790
Power: 23,0 kW
Voltage: 400 3/N V
Weight: 780 kg

High Temperature Chamber Furnaces HKT 160/18

Max. temperature: 1800°C
Inside dimensions (mm): W 500 x D 550 x H 550
Volume: 150 l
Outside dimensions (mm): W 1200 x D 975 x H 1940
Power: 21,0 kW
Voltage: 400 3/N V
Weight: 670 kg

High Temperature Chamber Furnaces HKT 220/18

Max. temperature: 1800°C
Inside dimensions (mm): W 500 x D 550 x H 800
Volume: 220 l
Outside dimensions (mm): W 1200 x D 975 x H 2190
Power: 27,0 kW
Voltage: 400 3/N V
Weight: 1000 kg

High Temperature Chamber Furnaces HKT 300/18

Max. temperature: 1800°C
Inside dimensions (mm): W 500 x D 1100 x H 550
Volume: 300 l
Outside dimensions (mm): W 1260 x D 1605 x H 1940
Power: 32,0 kW
Voltage: 400 3/N V
Weight: 1270 kg

High Temperature Chamber Furnaces HKT 500/18

Max. temperature:	1800°C
Inside dimensions (mm):	W 500 x D 1200 x H 800
Volume:	480 l
Outside dimensions (mm):	W 1260 x D 1705 x H 2190
Power:	58,0 kW
Voltage:	400 3/N V
Weight:	1310 kg