

## Technical Data High Temperature Combi Furnaces up to 1750

### High Temperature Combi Furnaces HTK 70/16 DB

<b>Max. temperature:</b>	<b>1600°C</b>
<b>Inside dimensions (mm):</b>	<b>W 400 x D 400 x H 400</b>
<b>Volume:</b>	<b>64 l</b>
<b>Outside dimensions (mm):</b>	<b>W 1100 x D 835 x H 1790</b>
<b>Power:</b>	<b>19,0 kW</b>
<b>Voltage:</b>	<b>400 3/N V</b>
<b>Weight:</b>	<b>550 kg</b>

### High Temperature Combi Furnaces HTK 160/16 DB

<b>Max. temperature:</b>	<b>1600°C</b>
<b>Inside dimensions (mm):</b>	<b>W 500 x D 550 x H 550</b>
<b>Volume:</b>	<b>150 l</b>
<b>Outside dimensions (mm):</b>	<b>W 1200 x D 975 x H 1940</b>
<b>Power:</b>	<b>25,0 kW</b>
<b>Voltage:</b>	<b>400 3/N V</b>
<b>Weight:</b>	<b>700 kg</b>

### High Temperature Combi Furnaces HTK 300/16 DB

<b>Max. temperature:</b>	<b>1600°C</b>
<b>Inside dimensions (mm):</b>	<b>W 500 x D 1100 x H 550</b>
<b>Volume:</b>	<b>300 l</b>
<b>Outside dimensions (mm):</b>	<b>W 1260 x D 1605 x H 1940</b>
<b>Power:</b>	<b>50,0 kW</b>
<b>Voltage:</b>	<b>400 3/N V</b>
<b>Weight:</b>	<b>1300 kg</b>

### High Temperature Combi Furnaces HTK 70/17 DB

**Max. temperature:** 1750°C  
**Inside dimensions (mm):** W 400 x D 400 x H 400  
**Volume:** 64 l  
**Outside dimensions (mm):** W 1100 x D 835 x H 1790  
**Power:** 19,0 kW  
**Voltage:** 400 3/N V  
**Weight:** 550 kg

### High Temperature Combi Furnaces HTK 160/17 DB

**Max. temperature:** 1750°C  
**Inside dimensions (mm):** W 500 x D 550 x H 550  
**Volume:** 150 l  
**Outside dimensions (mm):** W 1200 x D 975 x H 1940  
**Power:** 25,0 kW  
**Voltage:** 400 3/N V  
**Weight:** 700 kg

### High Temperature Combi Furnaces HTK 300/17 DB

**Max. temperature:** 1750°C  
**Inside dimensions (mm):** W 500 x D 1100 x H 550  
**Volume:** 300 l  
**Outside dimensions (mm):** W 1260 x D 1605 x H 1940  
**Power:** 50,0 kW  
**Voltage:** 400 3/N V  
**Weight:** 1300 kg